# Instructions for Use

## Intended use (Indications)

This Implant is a Class IIb single use device (non active surgical implant)

This Implant is intended to provide temporary internal fixation of bone.

## Information for use

The Patient should be made aware of this implants limitations to weight bearing and vigorous activity.

The Implant life is dependant on many surgical factors and should be determined by the surgeon.

## Contra indications

- Active or recent sepsis
- Atrophy or deformation of the bone
- Insufficient bone stock
- Infection
- Systemic and metabolic disorders
- Severe deformity
- Obesity
- Severe osteoporosis
- Voluntary or involuntary over activity or impact
- Tumour
- Allergy to implant materials
- Likelihood of falls
- Relevant disabilities
- Vascular insufficiency, muscular atrophy, neuromuscular disease, nervous disease in the region of the operative site
- Deficient or inadequate soft tissue in the operative site

## General warnings

- Do not use for any purpose outside the intended use of the device, as it may seriously affect the safety & function of the product.
- This Implant is not approved for attachment or fixation to the pedicles of the cervical, thoracic, or lumbar spine.
- Only for use by a qualified surgeon with prior knowledge of temporary internal fixation, who will determine the required procedure.
- Damage to, or modification of an implant, may cause stresses or defects that could result in implant failure.
- Under no circumstance must an implant be re-used. Each implant is designed for single use. Although an implant may appear undamaged, previous stresses to which it may have been subjected may cause failure if re-implanted.
- Any metallic device used in conjunction with this implant should be manufactured from a metal of a similar composition. Mixing materials may lead to metallic interactions such as corrosion, which will increase the risk of implant failure.
- Handle devices with care to prevent damage to surgical gloves.

- Patients who are overweight or who are physically very active will place relatively heavy demands on the implant. Such patients should be instructed on the limitations of the device to minimise the risk of premature failure.

## Specific product warnings

### Bone plates

- Unnecessary bending of the implant should be avoided as this may cause implant failure. If contouring of implants is necessary it should be performed without generating sharp bends, notching, or scratching.

### Bone screws

- Unnecessary tightening forces should be avoided as this may cause implant failure, and also make future removal of the implant difficult.

### Guide wire and cannulated screws combination

- Cannulations should be cleared during surgery to prevent the accumulation of debris and the potential binding and subsequent movement of the guide wire.
- It is recommended that continuous screening using image intensifiers is used during, guide wire insertion, and when cannulated devices are advanced over the guide wire. This will help eliminate any inadvertent guide wire advancement.

## Sterility

Implants manufactured by Global Orthopaedics (UK) Ltd, are supplied sterile or non-sterile as specified on the product labelling.

### Sterile information

- The integrity of packaging for each implant should be checked carefully to ensure that product sterility has not been compromised.
- Refer to package label for sterilization method and shelf life.
- Implants supplied sterile should not be re-sterilised.

### Non-sterile information

- It is the responsibility of the hospital to ensure that adequate procedures are taken to ensure product sterility prior to implantation.
- Materials used in Global Orthopaedics (UK) Ltd devices can be sterilised using steam sterilisation. Temperatures and times should be determined by the hospital.
- Warning. Hot items should not be handled or used as they can cause injury.
- Warning. If the packaging is damaged or open in any way prior to the sterilization process, or exposure to contamination is suspected, then this product should be disposed of. See Disposal below.
- The manufacturer and distributor take no responsibility for sterilisation or re-sterilisation of implants undertaken by the hospital.

## Disposal

- Dispose at end of life, or for any other reason, should be in accordance with hospital policy for Disposal of Contaminated Waste.

## Storage

- Store in a cool dry place between 0 and 30°C
- Keep packaging Dry
- Keep out of direct sunlight

## Post-operative care and patient information

- Patients who are overweight or who are physically very active will place relatively heavy demands on the implant. Such patients should be instructed on the limitations of the device to minimise the risk of premature failure.
- High stress and contact activities should be avoided during the healing process.
- The patient should present for post-operative review by the Surgeon at defined intervals.
- The patient should be encouraged to report any post-operative complications to the Surgeon.
- Post-operative wound infection can cause loosening of the implant, It is important that the surgeon and the care team provide the patient with information and guidance on how to maintain hygiene of the implant area.
- Smoking & non-steroidal anti-inflammatory drugs have both been shown to prolong the healing time of some fracture sites.

**Caution:** Please tell medical staff that you have an implant before you have a magnetic resonance image (MRI) taken. The metal in this implant can affect the quality of the images taken.

## Possible adverse effects

- Bone-screw loosening.
- Bone-screw site infection.
- Bone necrosis due to incorrect drilling or application, leading to early bone-screw loosening.
- Bone-screw bending or breakage due to over-stressing.
- Foreign body reaction to the implant caused by hypersensitivity against the metals. (Cobalt, Chrome, Titanium, Nickel). Can cause implant loosening or skin irritation.
- Improper reduction of the fracture fragments.
- Re-fracture through bone-screw holes.
- Pressure on the skin surrounding bone-screw sites due to improper dressing.
- Re-fracture due to removing the fixation too early.
- Non-union due to biological factors set at the time of injury or later.
- Problems due to non-compliance of patient.

## Symbol Definitions

SYMBOL	SYMBOL DEFINITION	SYMBOL	SYMBOL DEFINITION
	DO NOT REUSE / SINGLE USE ONLY		TEMPERATURE LIMITATION FOR STORAGE
	USE BY DATE YYYY/MM		KEEP AWAY FROM SUNLIGHT
	MANUFACTURERS BATCH CODE		KEEP DRY
	DATE OF MANUFACTURE		DO NOT RESTERILIZE
	STERILIZED USING IRRADIATION		PRODUCT NON STERILE
	CATALOGUE NUMBER / PRODUCT CODE		DO NOT USE IF PACKAGE IS DAMAGED
	CAUTION specific warnings associated with this product		CONSULT INSTRUCTIONS FOR USE
	MANUFACTURER	<b>SS</b>	Stainless steel
		<b>Ti</b>	Titanium

Please visit our web site to see the most up to date IFU

Other languages are available for download at our web site

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